



2015 Advanced Lithography :



Phase Measurements of EUV Mask Defects

Rene Claus



Yow-Gwo Wang, Markus Benk, Antoine Wojdyla, Alex Donoghue,
David Johnson, Kenneth Goldberg, Andrew Neureuther,
Patrick Naulleau, and Laura Waller

EUV Mask Defects

- Pits or particles in the substrate or multilayer will affect printing
- Can repair multilayer
- Can modify pattern
- Need to know what defect looks like
 - What is the phase/height?
 - What is the width?
 - Does it have amplitude?
- Want to measure using existing tools
 - AIMS tool




New Phase Retrieval Algorithm



Phase of EUV mask defect


Weak Object Assumption



- Consider a rough mirror (or mask)
 - Most of the light is reflected
 - Some of the light is scattered
- The electric field leaving the mask can be expressed as the sum of these components

$$E = 1 + E_s$$

$$I = |1 + E_s|^2 = 1 + 2\text{Re}\{E_s\} + \cancel{|E_s|^2}$$



- For most objects $\text{Scattering} \ll \text{DC}$
 - We can ignore *Scattering-Scattering*

Transfer Function

$$I = 1 + \text{Re}\{E_s\} * K_{re} + \text{Im}\{E_s\} * K_{im}$$

- For a weakly scattering object, we can express the intensity as two convolutions

$$F\{K_{re}\} = (P \cdot L) \star P + P \star (P \cdot L)$$

$$F\{K_{im}\} = (P \cdot L) \star P - P \star (P \cdot L)$$


P : pupil function, L : source shape

- $\text{Re}\{E_s\} \approx \text{amplitude}$, $\text{Im}\{E_s\} \approx \text{phase}$


$$E = (1 + A)e^{i\phi} \approx 1 + \underset{\substack{\uparrow \\ \text{Re}\{E_s\} \\ \text{(amplitude)}}}{A} + i\underset{\substack{\nwarrow \\ \text{Im}\{E_s\} \\ \text{(phase)}}}{\phi}$$

Coherent Transfer Function


Coherent Illumination




Partially Coherent Transfer Function



Inverting the Transfer Functions




Iterative Algorithm




- Ideally scattering-scattering term is zero
 - Subtract term \rightarrow becomes zero
- The iterative algorithm works reliably for defects even when they are large

SHARP



Zone Plate Microscope

- Pixel Size: 15nm
- NA: 0.33/4 (0.0825)
- Wavelength: 13.5nm
- Programmable illumination




Phase Contrast vs Standard Data

- Measured same defect with different pupils
- Qualitatively very different measurements
- If the recovered object is the same, it's not an artifact


$NA: 0.33/4$

$\sigma = 0.25$ illumination


Standard Pupil



Phase Contrast Pupil




Phase Contrast Imaging




- A phase contrast pupil has a phase shifting region
- Effectively switches phase and amplitude information

At Focus:

Standard




Phase Contrast



Y.G. Wang, "Enhancing defect detection with Zernike phase contrast in EUV multilayer blank inspection," SPIE Advanced Lithography (2015).



Standard Zone Plate Results




← $-4\mu\text{m}$ $0\mu\text{m}$ $4\mu\text{m}$ →

Phase (effective height): 2.5nm

Amplitude: 0.3





Phase Contrast Zone Plate Results



← $-4\mu m$ $0\mu m$ $4\mu m$ →


Phase (effective height): $2.5nm$

Amplitude: 0.3



Same Result

Comparing Results





Algorithm works for: complicated pupil and partial coherence ($\sigma = 0.25$)

Resolution Limited Defect


Very Small Native Defect

FWHM: 100nm , Height: 0.65nm , Absorption: 11%






Coherent results


Partial Coherence and Resolution



Source




Sub-Resolution Defect Printability




- Sub-resolution defects appear shorter and wider than they are under Coherent illumination
- Partially coherent measurements give a more accurate (higher resolution) phase result
- They print differently when patterned
- **Need to measure with partially coherently illumination or use higher NA inspection**

(Simulated Aerial Image)




16nm lines, dipole illumination
65nm defocus

Programmed Defect Study



Programmed Defects



- Can measure EUV defects smaller than what AFM can see
- Defects are not pure phase defects

Y.G. Wang, "Enhancing defect detection with Zernike phase contrast in EUV multilayer blank inspection," SPIE Advanced Lithography (2015).


T. Liang, "Growth and Printability of Multilayer Phase Defects on EUV Mask Blanks," JVSTB (2007).

Conclusion

- We've developed a new algorithm to get the phase and amplitude from measurements
 - Can use partially coherent illumination
 - Verified using phase contrast zone plate
- Small defects may be larger on the mask than seen with AIMS
 - Can use partially coherent illumination to get higher fidelity measurement of the mask phase
- Mask defects are not pure phase defects
 - Large defects have higher absorption

Acknowledgement






Phase of EUV mask defect

Thank you for your attention!

reneclaus@gmail.com

Berkeley
UNIVERSITY OF CALIFORNIA



Phase of EUV mask defect

Questions?

reneclaus@gmail.com

Berkeley
UNIVERSITY OF CALIFORNIA